


Product datasheet

Anti-ACBD4 antibody ab87176

1 图像

概述

产品名称	Anti-ACBD4抗体
描述	兔多克隆抗体to ACBD4
宿主	Rabbit
特异性	Recognises isoform 2
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Cow, Cat, Dog, Pig, Zebrafish 
免疫原	Synthetic peptide within residues: GTEKESPEPD CQKQFQAAVS VIQNLPKNGS YRPSYEMLR FYSYKQATM, corresponding to N terminal amino acids 2-51 of Human ACBD4 (NP_078998) Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Transfected 293T cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab87176** in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
----	-------	----


WB Use a concentration of 1 µg/ml. Detects a band of approximately 35 kDa (predicted molecular weight: 35 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

序列相似性

Contains 1 ACB (acyl-CoA-binding) domain.

图片



Western blot - ACBD4 antibody (ab87176)

Anti-ACBD4 antibody (ab87176) at 1 µg/ml +
Transfected 293T cell lysate at 10 µg

Secondary

anti-Rabbit IgG HRP at 1/50000 dilution

Predicted band size: 35 kDa

Observed band size: 35 kDa

Gel concentration 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors