


Product datasheet

Anti-ACBD4 antibody ab87176

1 图像

概述

<b>产品名称</b>	Anti-ACBD4抗体
<b>描述</b>	兔多克隆抗体to ACBD4
<b>特异性</b>	Recognises isoform 2
<b>经测试应用</b>	<b>适用于:</b> WB
<b>种属反应性</b>	<b>与反应:</b> Human <b>预测可用于:</b> Mouse, Rat, Rabbit, Horse, Cow, Cat, Dog, Pig, Zebrafish 
<b>免疫原</b>	Synthetic peptide within residues: GTEKESPEPD CQKQFQAAVS VIQNLPKNGS YRPSYEMLR FYSYKQATM, corresponding to N terminal amino acids 2-51 of Human ACBD4 (NP_078998) <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
<b>阳性对照</b>	Transfected 293T cell lysate.

性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
<b>存储溶液</b>	Preservative: None Constituents: 2% Sucrose, PBS
<b>纯度</b>	Immunogen affinity purified
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab87176** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
----	-------	----


WB Use a concentration of 1 µg/ml. Detects a band of approximately 35 kDa (predicted molecular weight: 35 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

### 序列相似性

Contains 1 ACB (acyl-CoA-binding) domain.

## 图片



Anti-ACBD4 antibody (ab87176) at 1 µg/ml +  
Transfected 293T cell lysate at 10 µg

### Secondary

anti-Rabbit IgG HRP at 1/50000 dilution

**Predicted band size** : 35 kDa

**Observed band size** : 35 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors