

## Product datasheet

# Anti-CIB3 antibody ab87784

### 1 图像

#### 概述

---

产品名称	Anti-CIB3抗体
描述	兔多克隆抗体to CIB3
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Sheep, Rabbit, Chicken, Guinea pig, Cow, Cat, Dog, Pig, Zebrafish



免疫原	Synthetic peptide corresponding to a region within N-terminal amino acids 36-85 (DLAPQLVPLD YTTCPDVKVP YELIGSMPEL KDNPFQRQRIA QVFSEDGDGH) of Human CIB3 (NP_473454).
-----	--

 Run BLAST with   Run BLAST with

阳性对照	Human fetal spleen lysate
------	---------------------------

#### 性能

---

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

---

Our [Abpromise guarantee](#) covers the use of **ab87784** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 22 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

## 相关性

The exact function of CIB3 is unknown, however, the gene product shares a high degree of sequence similarity with DNA-dependent protein kinase catalytic subunit-interacting protein 2 in human and mouse, and like them may bind the catalytic subunit of DNA-dependent protein kinases.

## 图片



Anti-CIB3 antibody (ab87784) at 1 µg/ml +  
Human fetal spleen lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000  
dilution

**Predicted band size: 22 kDa**

Gel concentration: 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors