

## Product datasheet

# Anti-RNF35 antibody ab88109

### 3 图像

#### 概述

产品名称	Anti-RNF35抗体
描述	小鼠多克隆抗体to RNF35
宿主	Mouse
经测试应用	适用于: WB, ICC/IF
种属反应性	与反应: Human
免疫原	Full length Human RNF35 protein (AAH60785.1)
阳性对照	K562 cell lysate, RNF35 transfected 293T cell lysate and HeLa cells

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab88109** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 29 kDa.
ICC/IF		Use a concentration of 10 µg/ml.

## 靶标

## 相关性


RNF35 belongs to the TRIM/RBCC family.

## 图片



Anti-RNF35 antibody (ab88109) at 1  $\mu\text{g/ml}$  +  
K562 cell lysate at 50  $\mu\text{g}$

**Predicted band size:** 29 kDa




**All lanes :** Anti-RNF35 antibody (ab88109) at  
1  $\mu\text{g/ml}$

**Lane 1 :** RNF35 transfected 293T cell lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25  $\mu\text{g}$  per lane.

**Predicted band size:** 29 kDa



ab88109 at 10 $\mu\text{g/ml}$  staining RNF35 in HeLa  
cells by Immunofluorescence.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Promise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors