


Product datasheet

Anti-NXF2 antibody ab88329

2 图像

概述

产品名称	Anti-NXF2抗体
描述	小鼠单克隆抗体to NXF2
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human, Recombinant fragment 预测可用于: Chimpanzee, Monkey 
免疫原	Recombinant fragment: DEIRITWRN RKPPERKMSQ NTQDGYTRNW FKVTIPYGIK YDKAWLMNSIQSHCSRFTP VDFHYVRNRA CFFV, corresponding to amino acids 95-169 of Human NXF2 with proprietary tag. Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	HeLa cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab88329** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 72 kDa.


应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.

靶标

相关性 NXF2 is involved in the export of mRNA from the nucleus to the cytoplasm.

细胞定位 Nuclear. Cytoplasmic.

图片




Western blot - NXF2 antibody (ab88329)

Anti-NXF2 antibody (ab88329) at 5 µg/ml +
HeLa cell lysate at 50 µg

Predicted band size : 72 kDa

Observed band size : 60 kDa

Additional bands at : 24 kDa. We are unsure as to the identity of these extra bands.



Western blot - NXF2 antibody (ab88329)

Anti-NXF2 antibody (ab88329) at 5 µg/ml +
recombinant protein fragment at 0.2 µg

Predicted band size : 72 kDa

Observed band size : 34 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors