

## Product datasheet

# Anti-NXF2 antibody ab88329

### 2 图像

#### 概述

产品名称	Anti-NXF2抗体
描述	小鼠单克隆抗体to NXF2
宿主	Mouse
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human, Recombinant fragment 预测可用于: Chimpanzee, Monkey
免疫原	Recombinant fragment: DEIRITWRN RKPPERKMSQ NTQDGYTRNW FKVTIPYGIK YDKAWLMNSIQSHCSRFTP VDFHYVRNRA CFFV, corresponding to amino acids 95-169 of Human NXF2 with proprietary tag. <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	HeLa cell lysate

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab88329** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 72 kDa.
ELISA		Use at an assay dependent dilution.

## 靶标


<b>相关性</b>	NXF2 is involved in the export of mRNA from the nucleus to the cytoplasm.
<b>细胞定位</b>	Nuclear. Cytoplasmic.

## 图片



Anti-NXF2 antibody (ab88329) at 5 µg/ml +  
HeLa cell lysate at 50 µg

**Predicted band size:** 72 kDa  
**Observed band size:** 60 kDa  
**Additional bands at:** 24 kDa. We are unsure  
as to the identity of these extra bands.



Anti-NXF2 antibody (ab88329) at 5 µg/ml +  
recombinant protein fragment at 0.2 µg

**Predicted band size:** 72 kDa  
**Observed band size:** 34 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors