

Product datasheet

Anti-FAM3C antibody ab88337

2 图像

概述

产品名称	Anti-FAM3C抗体
描述	小鼠单克隆抗体to FAM3C
宿主	Mouse
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human
免疫原	Recombinant fragment corresponding to amino acids 25 - 228 of Human FAM3C (AAH46932) with a proprietary tag
阳性对照	A431 cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: PBS, pH 7.2
纯度	Protein A purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab88337** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 25 kDa.
ELISA		Use at an assay dependent dilution.

靶标

功能	May be involved in retinal laminar formation. Promotes epithelial to mesenchymal transition.
组织特异性	Present in most secretory epithelia (at protein level).
序列相似性	Belongs to the FAM3 family.
细胞定位	Secreted. Cytoplasmic vesicle. Cytoplasmic in some cancer cells.

图片



Anti-FAM3C antibody (ab88337) at 1 µg/ml + immunogen at 0.2 µg


Secondary

Goat anti-Mouse IgG at 1/5000 dilution

Predicted band size: 25 kDa

Observed band size: 48 kDa

The observed band size may not correspond to the predicted protein molecular weight as the immunogen (a recombinant tagged protein) was used for the test lane.



Anti-FAM3C antibody (ab88337) at 1 µg/ml + A431 cell lysate at 25 µg

Secondary

Goat anti-Mouse IgG at 1/5000 dilution

Predicted band size: 25 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors