

## Product datasheet

# Anti-FAM3C antibody ab88337

### 2 图像

#### 概述

产品名称	Anti-FAM3C抗体
描述	小鼠单克隆抗体to FAM3C
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human
免疫原	Recombinant fragment corresponding to amino acids 25 - 228 of Human FAM3C (AAH46932) with a proprietary tag
阳性对照	A431 cell lysate

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: PBS, pH 7.2
纯度	Protein A purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab88337** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 25 kDa.
ELISA		Use at an assay dependent dilution.

## 靶标

<b>功能</b>	May be involved in retinal laminar formation. Promotes epithelial to mesenchymal transition.
<b>组织特异性</b>	Present in most secretory epithelia (at protein level).
<b>序列相似性</b>	Belongs to the FAM3 family.
<b>细胞定位</b>	Secreted. Cytoplasmic vesicle. Cytoplasmic in some cancer cells.

## 图片




Anti-FAM3C antibody (ab88337) at 1 µg/ml + immunogen at 0.2 µg

### Secondary

Goat anti-Mouse IgG at 1/5000 dilution

**Predicted band size** : 25 kDa

**Observed band size** : 48 kDa



Anti-FAM3C antibody (ab88337) at 1 µg/ml + A431 cell lysate at 25 µg

### Secondary

Goat anti-Mouse IgG at 1/5000 dilution

**Predicted band size** : 25 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors