

Product datasheet

Anti-Homez antibody ab88511

3 图像

概述

产品名称	Anti-Homez抗体
描述	小鼠多克隆抗体to Homez
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	Recombinant full length Homez (Human)
阳性对照	Human placenta lysate and Human esophagus tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab88511** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

IHC-P

应用说明	IHC-P: Use at a concentration of 3 µg/ml. WB: Use at a concentration of 1 µg/ml. Predicted molecular weight: 58 kDa.
------	---

Not yet tested in other applications.


Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性 Homez and members of ZHX family of zinc finger homeodomain factors are delineated as a subset within the superfamily of homeobox-containing proteins.

细胞定位 Nuclear


图片



Anti-Homez antibody (ab88511) at 1 $\mu\text{g/ml}$ +
Human placenta lysate at 50 μg

Predicted band size : 58 kDa

Observed band size : 58 kDa



All lanes : Anti-Homez antibody (ab88511) at
1 $\mu\text{g/ml}$

Lane 1 : transfected 293T cell lysate


Lane 2 : untransfected 293T cell lysate

Lysates/proteins at 25 μg per lane.

Predicted band size : 58 kDa

Observed band size : 58 kDa

Additional bands at : 40 kDa. We are
unsure as to the identity of these extra bands.



Purified ab88511 at 3ug/ml staining Homez in formalin-fixed, paraffin-embedded Human esophagus tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Homez antibody (ab88511)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors