

Product datasheet

Anti-CSH2 antibody ab88652

2 图像

概述

产品名称	Anti-CSH2抗体
描述	小鼠多克隆抗体to CSH2
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length human CSH2 protein (AAH22044).
阳性对照	Human placenta tissue lysate and CSH2 transfected 293T cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab88652** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 25 kDa.

靶标

功能 Produced only during pregnancy and is involved in stimulating lactation, fetal growth and

metabolism. Does not interact with GHR but only activates PRLR through zinc-induced dimerization.


序列相似性

Belongs to the somatotropin/prolactin family.

细胞定位

Secreted.

图片




Anti-CSH2 antibody (ab88652) at 1 $\mu\text{g/ml}$ +
human placenta tissue lysate at 50 μg

Developed using the ECL technique.

Predicted band size: 25 kDa

Observed band size: 23 kDa



All lanes : Anti-CSH2 antibody (ab88652) at
1 $\mu\text{g/ml}$

Lane 1 : cell lysate from CSH2 transfected
293T cells

Lane 2 : cell lysate from non transfected 293T
cells

Lysates/proteins at 25 μg per lane.

Developed using the ECL technique.

Predicted band size: 25 kDa

Observed band size: 22 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors