

Product datasheet

Anti-CPXCR1 antibody ab89527

3 图像

概述

产品名称	Anti-CPXCR1抗体
描述	小鼠多克隆抗体to CPXCR1
宿主	Mouse
经测试应用	适用于: WB, ICC/IF
种属反应性	与反应: Human
免疫原	Full length Human protein
阳性对照	WB: liver tissue lysate; ICC/IF: HeLa cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab89527** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 35 kDa.
ICC/IF		Use a concentration of 10 µg/ml.

靶标

组织特异性


Expressed in a variety of fetal tissues.

图片



Anti-CPXCR1 antibody (ab89527) at 1 μ g/ml
+ Human liver tissue lysate at 50 μ g

Predicted band size: 35 kDa




All lanes : Anti-CPXCR1 antibody (ab89527)
at 1 μ g/ml

Lane 1 : CPXCR1 transfected 293T cell
lysate

Lane 2 : Non-transfected lysate

Lysates/proteins at 25 μ g per lane.

Predicted band size: 35 kDa



ab89527 at 10 μ g/ml staining CPXCR1 in
PFA-fixed HeLa cells by immunofluorescence.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors