

Product datasheet

Anti-KRR1 antibody ab89601

3 图像

概述

产品名称	Anti-KRR1抗体
描述	小鼠多克隆抗体to KRR1
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length Human KRR1 (AAH26107.1)
阳性对照	Human pancreas tissue lysate. HeLa cell lysate. KRR1 transfected 293T cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab89601** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 44 kDa.

靶标

<b>功能</b>	Required for 40S ribosome biogenesis. Involved in nucleolar processing of pre-18S ribosomal RNA and ribosome assembly.
<b>序列相似性</b>	Belongs to the KRR1 family. Contains 1 KH domain.
<b>细胞定位</b>	Nucleus > nucleolus.


## 图片



Anti-KRR1 antibody (ab89601) at 1  $\mu$ g/ml +  
Human pancreas tissue lysate at 50  $\mu$ g

**Predicted band size:** 44 kDa


**Observed band size:** 44 kDa



Anti-KRR1 antibody (ab89601) at 1  $\mu$ g/ml +  
HeLa cell lysate at 50  $\mu$ g

**Predicted band size:** 44 kDa

**Observed band size:** 50 kDa



**All lanes :** Anti-KRR1 antibody (ab89601) at  
1  $\mu$ g/ml

**Lane 1 :** KRR1 transfected 293T cell lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25  $\mu$ g per lane.

**Predicted band size:** 44 kDa

**Observed band size:** 48 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors