

Product datasheet

Anti-Eph receptor B3 antibody ab95224

2 图像

概述

产品名称	Anti-Eph receptor B3抗体
描述	兔多克隆抗体to Eph receptor B3
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Human
免疫原	Recombinant Human Eph receptor B3 with His tag (NP_004434).
阳性对照	C6 cell line lysate and Mouse stomach tissue lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Ammonium Sulphate Precipitation
纯化说明	ab95224 is purified by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab95224** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/50 - 1/100. Predicted molecular weight: 110 kDa.

图标


<b>功能</b>	Receptor for members of the ephrin-B family. Binds to ephrin-B1 and -B2.
<b>组织特异性</b>	Ubiquitous.
<b>序列相似性</b>	Belongs to the protein kinase superfamily. Tyr protein kinase family. Ephrin receptor subfamily. Contains 2 fibronectin type-III domains. Contains 1 protein kinase domain. Contains 1 SAM (sterile alpha motif) domain.
<b>细胞定位</b>	Membrane.

## 图片



Anti-Eph receptor B3 antibody (ab95224) at  
1/50 dilution + C6 cell line lysate at 35 µg

**Predicted band size:** 110 kDa



Anti-Eph receptor B3 antibody (ab95224) at  
1/50 dilution + mouse stomach tissue lysate at  
35 µg

**Predicted band size:** 110 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors