

Product datasheet

Anti-ALPPL2 antibody ab99177

1 图像

概述

产品名称	Anti-ALPPL2抗体
描述	兔多克隆抗体to ALPPL2
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Chicken, Cow, Dog, Zebrafish 
免疫原	Synthetic peptide corresponding to a region within internal amino acids 431-480 ( SESGSPEYRQQSAVPLDGETHAGEDVAVFARGPQAHLVHGVQEQTFFIAHV ) of human ALPPL2 (NP_112603). <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
阳性对照	OVCAR-3 cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用


Our [Abpromise guarantee](#) covers the use of **ab99177** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 57 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

组织特异性	Trace amounts in the testis and thymus, and in elevated amounts in germ cell tumors.
序列相似性	Belongs to the alkaline phosphatase family.
细胞定位	Cell membrane.

## Anti-ALPPL2 antibody 图像



Anti-ALPPL2 antibody (ab99177) at 1  $\mu$ g/ml +  
OVCAR-3 cell lysate at 10  $\mu$ g

### Secondary

Goat Anti-Rabbit IgG Antibody (HRP-  
Conjugated)

**Predicted band size : 57 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors