

Product datasheet

Mouse monoclonal 2A 8F4 Anti-Rat IgG2a heavy chain (Biotin)  
ab99668

1 图像

概述

产品名称	小鼠单克隆抗体2A 8F4抗大鼠IgG2a heavy chain (Biotin)
宿主	Mouse
靶标种属	Rat
特异性	ab99668 reacts with the $\gamma$ 2a heavy chain (Fc) of Rat IgG2a.
经测试应用	适用于: ELISA, ICC/IF, Flow Cyt
偶联物	Biotin

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	2A 8F4
同种型	IgG1
轻链类型	kappa

应用


Our [Abpromise guarantee](#) covers the use of **ab99668** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.

应用	Ab评论	说明
ICC/IF		Use at an assay dependent concentration.
Flow Cyt		Use 1µg for 10 <sup>6</sup> cells. <a href="#">ab18434</a> -Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

## 图片



Top: (Negative control) ab99668 being used against BALB/c spleenocytes at 1 µg/10<sup>6</sup> cells with Streptavidin FITC.

Bottom: ab99668 staining BALB/c spleenocytes incubated with rat anti-mouse CD8 (IgG2a isotype) at 1 µg/10<sup>6</sup> cells with Streptavidin FITC.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors