

Product datasheet

Mouse monoclonal [31-7-4] Anti-Human IgG2 Fc ab99777

概述

| | |
|-------|---|
| 产品名称 | 小鼠单克隆抗体[31-7-4]抗人IgG2 Fc |
| 宿主 | Mouse |
| 靶标种属 | Human |
| 特异性 | Reacts with the Fc portion of the heavy chain of human IgG2 |
| 经测试应用 | 适用于: Dot blot, ELISA, WB, ICC/IF, IHC-Fr |
| 免疫原 | Other Immunogen Type corresponding to IgG2. |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| 存储溶液 | pH: 8.20 Constituent: 100% Borate buffered saline No preservatives or amine-containing buffer salts added. |
| 纯度 | Protein G purified |
| 克隆 | 单克隆 |
| 克隆编号 | 31-7-4 |
| 同种型 | IgG1 |
| 轻链类型 | kappa |

应用

Our [Abpromise guarantee](#) covers the use of **ab99777** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----------|------|-------------------------------------|
| Dot blot | | Use at an assay dependent dilution. |
| ELISA | | Use at an assay dependent dilution. |

| 应用 | Ab评论 | 说明 |
|--------|------|-------------------------------------|
| WB | | Use at an assay dependent dilution. |
| ICC/IF | | Use 1µg for 10 ⁶ cells. |
| IHC-Fr | | Use at an assay dependent dilution. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors