

Product datasheet

Mouse monoclonal HP6014 Anti-Human IgG2 Fd gamma (Biotin) ab99791

概述

产品名称	小鼠单克隆抗体HP6014抗人IgG2 Fd gamma (Biotin)
宿主	Mouse
靶标种属	Human
特异性	ab99791 reacts with the Fd portion of the heavy chain of Human IgG2 as demonstrated by ELISA
经测试应用	适用于: ICC/IF, ELISA, WB, Dot blot, ICC
偶联物	Biotin

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS
纯度	Protein G purified
纯化说明	Purified from ascites fluid.
克隆	单克隆
克隆编号	HP6014
同种型	IgG1
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab99791** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		Use 1µg for 10 ⁶ cells.

应用	Ab评论	说明
ELISA		1/1000 - 1/2000.
WB		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution.
ICC		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors