

Product datasheet

Mouse monoclonal [A9604D2] Anti-Human IgA2 Fc (Phycoerythrin) ab99803

概述

产品名称	小鼠单克隆抗体[A9604D2]抗人IgA2 Fc (Phycoerythrin)
宿主	Mouse
靶标种属	Human
特异性	ab99803 reacts with the Fc portion of the heavy chain of IgA ₂ .
经测试应用	适用于: ELISA, Dot blot, WB, ICC/IF, IHC-Fr
免疫原	Other Immunogen Type corresponding to IgA ₂ .
偶联物	Phycoerythrin. Ex: 488nm, Em: 575nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS, Sucrose
纯度	Protein G purified
克隆	单克隆
克隆编号	A9604D2
同种型	IgG1
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab99803** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution. for Dot and slot immunoblotting.

应用	Ab评论	说明
WB		Use at an assay dependent dilution. Predicted molecular weight: 37 kDa.
ICC/IF		Use 0.2µg for 10 ⁶ cells.
IHC-Fr		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors