

## Product datasheet

# Mouse monoclonal [HP6061] Anti-Human IgE Fc ab99833

### 概述

产品名称	小鼠单克隆抗体[HP6061]抗人IgE Fc
宿主	Mouse
靶标种属	Human
特异性	ab99833 reacts with the Fc portion of the heavy chain of Human IgE as demonstrated by ELISA.
经测试应用	<b>适用于: ELISA</b>
免疫原	Other Immunogen Type corresponding to IgE.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 8.20 Constituent: 100% Borate buffered saline  No preservatives or amine-containing buffer salts added.
纯度	Purified IgM
克隆	单克隆
克隆编号	HP6061
同种型	IgM

### 应用

Our [Abpromise guarantee](#) covers the use of **ab99833** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use a concentration of 2.5 - 10 µg/ml. This antibody is best used as a capture antibody.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors