

Product datasheet

Mouse monoclonal [HP6029] Anti-Human IgE Fc (HRP) ab99836

1 References

概述

产品名称	小鼠单克隆抗体[HP6029]抗人IgE Fc (HRP)
宿主	Mouse
靶标种属	Human
特异性	ab99836 reacts with the Fc portion of the heavy chain of Human IgE as demonstrated by ELISA
经测试应用	适用于: Dot blot, ELISA, IHC-Fr, WB, ICC
免疫原	Other Immunogen Type corresponding to IgE.
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: None Constituents: 50% PBS, 50% Glycerol, pH 7.4
纯度	Protein A purified
纯化说明	Source: Ascites fluid
克隆	单克隆
克隆编号	HP6029
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab99836** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Dot blot		Use at an assay dependent dilution.
ELISA		1/1000 - 1/2000.

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution.
ICC		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors