

Product datasheet

Mouse monoclonal [SB139e] Anti-Armenian Hamster IgG2+IgG3 heavy chain (Biotin) ab99842

1 References

概述

产品名称	小鼠单克隆抗体[SB139e]抗亚美尼亚仓鼠IgG2+IgG3 heavy chain (Biotin)
宿主	Mouse
靶标种属	Armenian hamster
特异性	Reacts with the $\gamma 2$ heavy chain of Syrian and Armenian Hamster IgG2 as well as with the $\gamma 3$ heavy chain of Armenian Hamster IgG3.
经测试应用	适用于: ELISA
免疫原	Other Immunogen Type corresponding to IgG2 + IgG3.
偶联物	Biotin

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS
纯度	IgG fraction
克隆	单克隆
克隆编号	SB139e
同种型	IgG1
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab99842** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		1/4000 - 1/8000.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors