

Anti-alpha MSH antibody ab50013

概述

产品名称	Anti-alpha MSH抗体
描述	兔多克隆抗体to alpha MSH
特异性	Anti-a-MSH cross-reacts with des-Ac-a-MSH and [Val-OH13]-a-MSH. Minimal cross-reactivity is observed with human ACTH (1-39) and ACTH (1-24).
经测试应用	适用于: RIA
种属反应性	与反应: Human
免疫原	Synthetic human MSH alpha peptide conjugated to bovine serum albumin. The mid-region sequence of MSH alpha (amino acids 4-10) is identical to the ACTH and MSH beta sequences and highly conserved in the MSH gamma sequence.
常规说明	Affinity constant K_a : at least 1×10^9 L/M

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 15mM Sodium Azide Constituents: 0.01M PBS, pH 7.4
纯度	Ion Exchange Chromatography
纯化说明	Whole antiserum is fractionated and then further purified by ion-exchange chromatography to provide the IgG fraction of antiserum that is essentially free of other rabbit serum proteins.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab50013** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

应用	Ab评论	说明
RIA		
应用说明	RIA: 1/3000 (2nd antibody and polyethylene glycol method) using 5-7pg of ¹²⁵ I-a-MSH. Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user.	
靶标		
相关性	Melanocyte-stimulating hormone belongs to a group called the melanocortins. This group includes ACTH, alpha-MSH, beta-MSH and gamma-MSH; these peptides are all cleavage products of a large precursor peptide called pro-opiomelanocortin (POMC). Alpha-MSH is the most important melanocortin for pigmentation. It increases the pigmentation of skin by increasing melanin production in melanocytes.	
细胞定位	Secreted	

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors