

Product datasheet

Anti-beta Defensin 1 antibody ab115813

1 References 1 图像

概述

产品名称	Anti-beta Defensin 1抗体
描述	兔多克隆抗体to beta Defensin 1
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	Synthetic peptide corresponding to a sequence at the C terminal of Human beta Defensin 1.
阳性对照	Human lung cancer tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservatives: 0.025% Thimerosal (merthiolate), 0.025% Sodium azide Constituents: 2.5% BSA, 0.45% Sodium chloride, 0.1% Dibasic monohydrogen sodium phosphate
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab115813** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
----	-------	----

WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 7 kDa. The detection limit for beta Defensin 1 is approximately 1ng/lane under non-reducing and reducing conditions.
----	--	--


应用	Ab 评论	说明
----	-------	----

IHC-P Use a concentration of 1 µg/ml. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

## 靶标

功能	Has bactericidal activity.
组织特异性	Plasma.
序列相似性	Belongs to the beta-defensin family.
细胞定位	Secreted.

## 图片



ab115813 at 1µg/ml staining beta Defensin 1 in by Paraffin-embedded Human Lung cancer tissue by Immunohistochemistry,

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-beta Defensin 1 antibody (ab115813)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors