

Product datasheet

Anti-beta Defensin 1 antibody ab14421

2 References

概述

产品名称	Anti-beta Defensin 1抗体
描述	兔多克隆抗体to beta Defensin 1
宿主	Rabbit
经测试应用	适用于: ELISA, WB, IHC-Fr
种属反应性	与反应: Human
免疫原	Synthetic peptide: YRGKAKC(ACM)G , corresponding to amino acids 28-34 of Human beta 1 Defensin.
	Run BLAST with Run BLAST with

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: Whole Serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab14421** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration. Predicted molecular weight: 8 kDa.

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent concentration. PubMed: 21056450

靶标

功能	Has bactericidal activity.
组织特异性	Plasma.
序列相似性	Belongs to the beta-defensin family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors