

Product datasheet

Anti-intestinal alkaline phosphatase antibody [AAP.1] ab18111

2 References

概述

产品名称	Anti-intestinal alkaline phosphatase抗体[AAP.1]
描述	小鼠单克隆抗体[AAP.1] to intestinal alkaline phosphatase
宿主	Mouse
经测试应用	适用于: ELISA
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human intestinal alkaline phosphatase. D98/AH-2(HeLa) cells, which produce alkaline phosphatase ectopically.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: None Constituents: PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	AAP.1
同种型	IgG2a

应用

Our [Abpromise guarantee](#) covers the use of **ab18111** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.

靶标

序列相似性	Belongs to the alkaline phosphatase family.
细胞定位	Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors