

## Product datasheet

# Anti-intestinal alkaline phosphatase antibody [AAP.1] ab18111

## 2 References

### 概述

产品名称	Anti-intestinal alkaline phosphatase抗体[AAP.1]
描述	小鼠单克隆抗体[AAP.1] to intestinal alkaline phosphatase
经测试应用	适用于: ELISA
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human intestinal alkaline phosphatase. D98/AH-2(HeLa) cells, which produce alkaline phosphatase ectopically.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: None Constituents: PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	AAP.1
同种型	IgG2a

### 应用

Our [Abpromise guarantee](#) covers the use of **ab18111** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.

### 靶标

序列相似性	Belongs to the alkaline phosphatase family.
细胞定位	Cell membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors