

Human Cytokine Antibody Array - Membrane (51 targets) ab169815

1 图像

概述

产品名称	人Cytokine抗体阵列-膜(51 targets)
样品类型	Cell culture supernatant, Saliva, Milk, Urine, Serum, Plasma, Cell culture extracts, Other biological fluids, Whole Blood, Tissue Extracts, Cell Lysate, Cell culture media
检测类型	Semi-quantitative
种属反应性	与反应: Human
产品概述	ab169815 is for the simultaneous detection of 51 Human Cytokines in 4 or 8 samples. Suitable for all liquid sample types.

Targets: Adipsin, BCAM, CD30, CD40, Fc gamma RIIB/C, Ferritin, FLRG, Follistatin, Furin, Galectin-7, GDF-15, Growth Hormone, IL-10 R alpha, IL-22, IL-28A, IL-29, IL-31, Insulin, Luteinizing hormone, LIMP2, LYVE-1, Marapsin, MICA, MICB, MMP-2, MMP-7, MMP-8, MMP-10, NCAM-1, Nidogen-1, NrCAM, NRG1-beta 1, Osteopontin, PAI-1, Platelet factor 4, PSA-total, RAGE, RANK, Resistin, SAA, Siglec-9, TACE, TIM-1, TRAIL R2, Trappin-2, TREM-1, TSH, TSLP, VCAM-1, VEGF-C, XEDAR

[Learn more about membrane antibody arrays](#)

说明

Abcam has not and does not intend to apply for the REACH Authorisation of customers' uses of products that contain European Authorisation list (Annex XIV) substances.

It is the responsibility of our customers to check the necessity of application of REACH Authorisation, and any other relevant authorisations, for their intended uses.

经测试应用

适用于: Multiplex Protein Detection

性能

存放说明

Store at -20°C. Please refer to protocols.

组件	1 x 4 Membranes	1 x 8 Membranes	1 x 2 Membranes
1,000X HRP-Conjugated Streptavidin	1 x 50µl	1 x 50µl	1 x 50µl
1X Blocking Buffer	1 x 25ml	2 x 25ml	1 x 25ml

组件	1 x 4 Membranes	1 x 8 Membranes	1 x 2 Membranes
20X Wash Buffer I	1 x 10ml	1 x 20ml	1 x 10ml
20X Wash Buffer II	1 x 10ml	1 x 20ml	1 x 10ml
2X Cell Lysis Buffer	1 x 10ml	1 x 16ml	1 x 10ml
8-Well Incubation Tray (with Lid)	1 unit	1 unit	1 unit
Biotinylated Antibody Cocktail (C9)	1 x 2 vials	1 x 4 vials	1 x 2 vials
Detection Buffer C	1 x 1.5ml	1 x 2.5ml	1 x 1.5ml
Detection Buffer D	1 x 1.5ml	1 x 2.5ml	1 x 1.5ml
Human Cytokine Antibody Array Membranes (C9)	1 x 4 units	1 x 8 units	1 x 2 units


应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab169815于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Multiplex Protein Detection		Use at an assay dependent concentration.

图片



Typical images obtained with Abcam Human Cytokine Antibody Array. These membranes were probed with conditioned media from two different cell lines. Membranes were exposed to film at room temperature for 1 minute.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors