

Phosphatase Inhibitor Cocktail I ab201112

4 References

概述

产品名称

Phosphatase抑制剂Cocktail I

产品概述

Abcam's phosphatase inhibitor cocktails inhibit a number of phosphatases to preserve the phosphorylation state of the proteins in your sample.

Target phosphatases - acid phosphatases, alkaline phosphatases, serine/threonine phosphatases, protein tyrosine phosphatases

Inhibitors included - sodium fluoride, sodium orthovanadate, sodium pyrophosphate decahydrate, β -glycerophosphate

Form: Supplied in water (1 ml)

Preparation: The recommended starting dilution is 1:49

Visit our [inhibitor cocktails page](#) for more information on our complete range of convenient and easy to use protease and phosphatase inhibitor cocktails.

说明

Individual components and concentrations:

Sodium fluoride: acid phosphatase and serine/threonine phosphatase inhibitor (250 mM)

Sodium orthovanadate: protein tyrosine phosphatase and alkaline phosphatase inhibitor (50 mM)

Sodium pyrophosphate decahydrate: serine/threonine phosphatase inhibitor (50 mM)

β -Glycerophosphate: serine/threonine phosphatase inhibitor (50 mM)

性能

存放说明

Store at -20°C. Please refer to protocols.

组件	1 ml
Phosphatase Inhibitor Cocktail I	1 x 1ml

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors