


### U0126, Selective MKK inhibitor ab120241

**20 References**   **2 图像**

#### 概述


产品名称	U0126, Selective MKK抑制剂
描述	Selective MKK抑制剂
生物学描述	Selective non-competitive inhibitor of MAP kinase kinase (MKK). Inhibits MKK1, MKK2 (IC <sub>50</sub> values of 0.07 and 0.06 μM, respectively) and MKK5 with little or no effect on a wide range of other kinases. Centrally active following systemic administration <i>in vivo</i> .
纯度	> 98%
CAS编号	109511-58-2
化学结构	

#### 性能

化学名称	1,4-Diamino-2,3-dicyano-1,4-bis(2-aminophenylthio)butadiene
分子量	380.49
分子式	C <sub>18</sub> H <sub>16</sub> N <sub>6</sub> S <sub>2</sub>
PubChem识别号	3006531
存放说明	Store at +4°C. Store under desiccating conditions. The product can be stored for up to 12 months.
溶解度概述	Soluble in DMSO to 100 mM
处理	<p>Wherever possible, you should prepare and use solutions on the same day. However, if you need to make up stock solutions in advance, we recommend that you store the solution as aliquots in tightly sealed vials at -20°C. Generally, these will be useable for up to one month. Before use, and prior to opening the vial we recommend that you allow your product to equilibrate to room temperature for at least 1 hour.</p> <p>Need more advice on solubility, usage and handling? Please visit our <a href="#">frequently asked questions (FAQ) page</a> for more details.</p>


**SMILES**
N/C(Sc1ccccc1N)=C(/C#N)C(C#N)=C(/N)Sc2ccccc2N
**来源**

Synthetic

**图片**

Chemical Structure - U0126, Selective MKK inhibitor  
(ab120241)

2D chemical structure image of ab120241, U0126, Selective MKK inhibitor



Western blot - U0126, Selective MKK inhibitor  
(ab120241)

Performed under reducing conditions; Samples loaded at 30 µg/lane. Membrane blocking and secondary antibody incubation performed in 5% milk/PBST. Primary antibody incubation performed in 5% BSA/PBST.

Lane 1: Marker

Lane 2: A431 lysate starved

Lane 3: A431 lysate starved treated with 100 ng/mL of EGF for 30 min

Lane 4: A431 lysate starved and pretreated for 2 hours with 10µM U0126 (ab120241)

Lane 5: A431 lysate starved and pretreated for 2 hours with 30µM Ly294002 (**ab120243**)

Primary: AKT/MAPK signaling pathway antibody cocktail at 1/250 dilution

Secondary: HRP conjugated Goat anti-Rabbit secondary antibody at 1:10,000 dilution.

Observed bands:

RSK1p90 (pS380) band size: 83 kDa

AKT1 (pS473) band size: 56 kDa

ERK1 (pY204) and ERK2 (pY187) band size: 43 kDa

RPS6 (pS235/236) band size: 31 kDa

Rab11a band size: 24 kDa

Unknown band: 130kDa

Performed under reducing conditions; Samples I

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES, NOT FOR USE IN HUMANS"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors
- Abcam biochemicals are novel compounds and we have not tested their biological activity in house. Please use the literature to identify how to use these products effectively. If you require further assistance please contact the scientific support team