

# Blocking Buffer 10X ab210904

## 2 References

### 概述

产品名称	Blocking缓冲液10X
产品概述	<p>Blocking Buffer 10X is recommended for use with the Matched Antibody Pair kits during ELISA development. The Blocking Buffer 10X is a BSA based buffer (10 %w/v) that is designed to reduce background interference due to non-specific binding.</p> <p>The 10X Blocking Buffer should be diluted to 1X using 1X PBS before use.</p> <p>See our <a href="#">website</a> for a protocol.</p> <p>The Blocking Buffer 10X is also available as part of our Matched Antibody Pair Kit Accessory Pack (<a href="#">ab210905</a>).</p>

### 性能

存放说明	Store at +4°C. Please refer to protocols.
------	---

组件	50 ml
Blocking Buffer 10X (ab210904)	1 x 50ml

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors