

Pepsin Solution (Antigen Retrieval) ab64201

1 图像

概述

产品名称	Pepsin Solution (Antigen Retrieval)
经测试应用	适用于: IHC-P
常规说明	<p>Pepsin for protease-induced antigen retrieval (PIER) during IHC procedures (ab64201). Ready-to-use formulation. Tris buffer pH 2.0</p> <p>Staining Protocol:</p> <p>Incubate tissue section for 10 minutes at 37°C.</p> <p>Rinse and proceed with immunostaining.</p> <p>Pepsin is a commonly used digestive enzyme for immunohistochemical procedures.</p> <p>Other products for IHC</p> <p>Other enzymes for antigen retrieval include: <u>Proteinase K Antigen Retrieval Solution ab64220</u> and <u>Trypsin Antigen Retrieval Kit ab259256</u>.</p> <p>Find more kits and reagents for antigen retrieval, blocking, signal amplification, visualization, counterstaining, and mounting in the <u>IHC kits and reagents guide</u>.</p>

性能

形式	Liquid
存放说明	Store at +4°C.
存储溶液	<p>pH: 1.90</p> <p>Preservative: 0.5% Sodium azide</p> <p>Constituents: 0.25% Pepsin, 0.22% Polyoxyethylene lauryl ether, 1.5% Hydrochloric acid</p>
Reagent说明	Pepsin is a commonly used digestive enzyme for immunohistochemical procedures.
功能	Shows particularly broad specificity; although bonds involving phenylalanine and leucine are preferred, many others are also cleaved to some extent.
序列相似性	Belongs to the peptidase A1 family.
细胞定位	Secreted.
细胞定位	Secreted.

应用


The Abpromise guarantee

Abpromise™ 承诺保证使用ab64201于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.

图片



ab64201 Pepsin Solution used with Laminin on Human Kidney

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Pepsin Solution (ab64201)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors