

Stop Solution 1X ab210900

概述

产品名称 Stop Solution 1X

产品概述 Stop Solution is recommended for use with the Matched antibody pair kits for sandwich ELISA. See our [website](#) for a protocol.

Stop Solution is supplied ready to use and is required for the stopping the action of horseradish peroxidase on the tetramethylbenzidine (TMB) substrate in an ELISA. When the TMB-peroxidase reaction is stopped using Stop Solution for TMB Substrates, the chromogen changes from blue to yellow and hence the absorbance from 620 nm to 450 nm.

The Stop Solution is also available as part of our Matched Antibody Pair Kit Accessory Pack ([ab210905](#)).

性能

存放说明 Store at +4°C. Please refer to protocols.

组件	120 ml
Stop Solution (ab210900)	1 x 120ml

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors