



Product datasheet

Universal antibody diluent ab79995

1 References

概述

产品名称	Universal抗体稀释剂
经测试应用	适用于: IHC-P
常规说明	Universal antibody diluent is a PBS formulation at pH 7.3 with protein carrier and preservative that has been developed for stability during long-term storage and for improved antibody performance. Protease and IgG-free.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.30 Preservative: 0.1% Sodium azide Constituents: PBS, Carrier protein  Proprietary modified PBS formulation; combination of buffers, staining enhancers, stabilizers and preservative

应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab79995于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors