

Annexin V/ANXA5 overexpression 293T lysate (whole cell)
ab76216

2 图像

概述

产品名称	Annexin V/ANXA5 overexpression 293T裂解物(whole cell)
常规说明	ab76216 is a 293T cell transfected lysate in which Human Annexin V/ANXA5 has been transiently over-expressed using a pCMV-Annexin V/ANXA5 plasmid. The lysate is provided in 1X Sample Buffer.
经测试应用	适用于: WB

性能


Mycoplasma free	Yes
形式	Liquid
存放说明	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 2% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.79% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Disease: Pregnancy loss, recurrent, 3 Domain: The [IL]-x-C-x-x-[DE] motif is a proposed target motif for cysteine S-nitrosylation mediated by the iNOS-S100A8/A9 transnitrosylase complex. A pair of annexin repeats may form one binding site for calcium and phospholipid. Function: This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is involved in the blood coagulation cascade. PTM: S-nitrosylation is induced by interferon-gamma and oxidatively-modified low-density lipoprotein (LDL(ox)) possibly implicating the iNOS-S100A8/9 transnitrosylase complex. Similarity: Belongs to the annexin family. Contains 4 annexin repeats.

应用


The Abpromise guarantee Abpromise™承诺保证使用ab76216于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent concentration.



ab76216 at 15µg/lane on an SDS-PAGE gel.



All lanes : Anti-Annexin V/ANXA5 antibody [1F4-1A5] (**ab54775**)
at 1/500 dilution

Lane 1 : Annexin V/ANXA5 overexpression 293T lysate (whole cell) (ab76216)

Lane 2 : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugated at
1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors