

C13orf31 overexpression 293T lysate (whole cell) ab215389

2 图像

概述

产品名称	C13orf31 overexpression 293T裂解物(whole cell)
常规说明	ab215389 is a 293T cell transfected lysate in which Human C13orf31 has been transiently over-expressed using a pCMV-C13orf31 plasmid. The lysate is provided in 1X Sample Buffer.
经测试应用	适用于: WB

性能

Mycoplasma free	Yes
形式	Liquid
存放说明	Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituents: 2% Sodium lauryl sulfate, 10% Glycerol (glycerin, glycerine), 0.01% Bromophenol blue, 0.79% Tris HCl, 2.34% Beta mercaptoethanol
背景	Similarity: Belongs to the UPF0124 family.


应用

The Abpromise guarantee Abpromise™承诺保证使用ab215389于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent concentration. Predicted molecular weight: 48 kDa.

图片



Western blot - C13orf31 overexpression 293T lysate (whole cell) (ab215389)


All lanes : Anti-C13orf31 antibody (**ab168238**)

Lane 1 : C13orf31 overexpression 293T lysate (whole cell) (ab215389)

Lane 2 : Untransfected 293T cell lysate

Lysates/proteins at 15 µl per lane.

Predicted band size: 48 kDa



SDS-PAGE - C13orf31 overexpression 293T lysate (whole cell) (ab215389)

ab215389 on an SDS Page gel.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors