

## Product datasheet

# Carboxypeptidase A overexpression 293T lysate (whole cell) ab94084

## 2 图像

### 概述

产品名称	Carboxypeptidase A overexpression 293T裂解物(whole cell)
常规说明	ab94084 is a 293T cell transfected lysate in which Human Carboxypeptidase A has been transiently over-expressed using a pCMV-Carboxypeptidase A plasmid. The lysate is provided in 1X Sample Buffer.
经测试应用	适用于: WB

### 性能

Mycoplasma free	Yes
形式	Liquid
存放说明	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: Carboxypeptidase that catalyzes the release of a C-terminal amino acid, but has little or no action with -Asp, -Glu, -Arg, -Lys or -Pro. Similarity: Belongs to the peptidase M14 family.

### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab94084于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		Use at an assay dependent dilution.

### 图片



SDS-PAGE - Carboxypeptidase A overexpression  
293T lysate (whole cell) (ab94084)

ab94084 at 15µg/lane on an SDS-PAGE gel.



Western blot - Carboxypeptidase A overexpression  
293T lysate (whole cell) (ab94084)

**All lanes** : Anti-Carboxypeptidase A antibody (**ab89592**) at 1/500 dilution

**Lane 1** : Carboxypeptidase A overexpression 293T lysate (whole cell) (ab94084)

**Lane 2** : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

#### Secondary

**All lanes** : Goat Anti-mouse IgG (H and L) HRP conjugated at 1/2500 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors