

CHRAC-17 overexpression 293T lysate (whole cell) ab94038

2 图像

概述

产品名称	CHRAC-17 overexpression 293T裂解物(whole cell)
常规说明	ab94038 is a 293T cell transfected lysate in which human CHRAC-17 has been transiently over-expressed using a pCMV-CHRAC-17 plasmid. The lysate is provided in 1X Sample Buffer. Protein previously known as DNA Polymerase epsilon.
经测试应用	适用于: WB

性能

Mycoplasma free	Yes
形式	Liquid
存放说明	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: Forms a complex with DNA polymerase epsilon subunit CHRAC1 and binds naked DNA, which is then incorporated into chromatin, aided by the nucleosome-remodeling activity of ISW/SNF2H and ACF1. Tissue specificity: Expressed in all tissues tested, including, heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas.


应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab94038于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		Use at an assay dependent dilution.

图片



SDS-PAGE - CHRAC-17 overexpression 293T lysate
(whole cell) (ab94038)

ab94038 at 15 µg/lane on an SDS-PAGE gel.



Western blot - CHRAC-17 overexpression 293T
lysate (whole cell) (ab94038)

All lanes : Anti-CHRAC-17 antibody ([ab88521](#)) at 1/500 dilution

Lane 1 : DNA Polymerase epsilon 293T Transfected Lysate -
Positive Control (ab94038)

Lane 2 : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugated at
1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors