

Cytochrome P450 2A6 overexpression 293T lysate (whole cell)  
ab94284

2 图像

概述

产品名称	Cytochrome P450 2A6 overexpression 293T裂解物(whole cell)
常规说明	ab94284 is a 293T cell transfected lysate in which Human Cytochrome P450 2A6 has been transiently over-expressed using a pCMV-Cytochrome P450 2A6 plasmid. The lysate is provided in 1X Sample Buffer.
经测试应用	适用于: WB

性能

Mycoplasma free	Yes
形式	Liquid
存放说明	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: Exhibits a high coumarin 7-hydroxylase activity. Can act in the hydroxylation of the anti-cancer drugs cyclophosphamide and ifosphamide. Competent in the metabolic activation of aflatoxin B1. Constitutes the major nicotine C-oxidase. Possesses low phenacetin O-deethylation activity. Tissue specificity: Liver. Similarity: Belongs to the cytochrome P450 family.


应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab94284于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent dilution.

图片



ab94284 at 15µg/lane on an SDS-PAGE gel.



**All lanes :** Anti-CYP2A6 antibody ([ab89633](#)) at 1/500 dilution

**Lane 1 :** Cytochrome P450 2A6 overexpression 293T lysate  
(whole cell) (ab94284)

**Lane 2 :** 293T Non Transfected Lysate

Lysates/proteins at 25 µg per lane.

#### Secondary

**All lanes :** Goat Anti-mouse IgG (H and L) HRP conjugated at  
1/2500 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors