

# Cytohesin 1 overexpression 293T lysate (whole cell) ab94252

## 2 图像

### 概述

产品名称	Cytohesin 1 overexpression 293T裂解物(whole cell)
常规说明	ab94252 is a 293T cell transfected lysate in which Human Cytohesin 1 has been transiently over-expressed using a pCMV-Cytohesin 1 plasmid. The lysate is provided in 1X Sample Buffer.
经测试应用	适用于: WB

### 性能

Mycoplasma free	Yes
形式	Liquid
存放说明	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: Promotes guanine-nucleotide exchange on ARF1 and ARF5. Promotes the activation of ARF through replacement of GDP with GTP. Tissue specificity: Ubiquitous. Similarity: Contains 1 PH domain. Contains 1 SEC7 domain.


### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab94252于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		Use at an assay dependent dilution.

### 图片



SDS-PAGE - Cytohesin 1 overexpression 293T lysate (whole cell) (ab94252)

ab94252 at 15µg/lane on an SDS-PAGE gel



Western blot - Cytohesin 1 overexpression 293T lysate (whole cell) (ab94252)

**All lanes** : Anti-Cytohesin 1 antibody ([ab56511](#)) at 1/500 dilution

**Lane 1** : Cytohesin 1 overexpression 293T lysate (whole cell) (ab94252)

**Lane 2** : 293T Non Transfected Lysate

Lysates/proteins at 25 µg per lane.

**Secondary**

**All lanes** : Goat Anti-mouse IgG (H and L) HRP conjugated at 1/2500 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors