

HDAC6 overexpression 293T lysate (whole cell) ab94173

2 图像

概述

产品名称	HDAC6 overexpression 293T裂解物(whole cell)
常规说明	ab94173 is a 293T cell transfected lysate in which Human HDAC6 has been transiently over-expressed using a pCMV-HDAC6 plasmid. The lysate is provided in 1X Sample Buffer.
经测试应用	适用于: WB

性能


Mycoplasma free	Yes
形式	Liquid
存放说明	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: Responsible for the deacetylation of lysine residues on the N-terminal part of the core histones (H2A, H2B, H3 and H4). Histone deacetylation gives a tag for epigenetic repression and plays an important role in transcriptional regulation, cell cycle progression and developmental events. Histone deacetylases act via the formation of large multiprotein complexes (By similarity). Plays a central role in microtubule-dependent cell motility via deacetylation of tubulin. Similarity: Belongs to the histone deacetylase family. HD type 2 subfamily. Contains 1 UBP-type zinc finger. PTM: Ubiquitinated. Its polyubiquitination however does not lead to its degradation. Sumoylated in vitro.

应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab94173于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent dilution.



SDS-PAGE - HDAC6 overexpression 293T lysate  
(whole cell) (ab94173)

ab94173 at 15µg/lane on an SDS-PAGE gel.



Western blot - HDAC6 overexpression 293T lysate  
(whole cell) (ab94173)

**All lanes :** Anti-HDAC6 antibody (**ab56926**) at 1/500 dilution

**Lane 1 :** HDAC6 overexpression 293T lysate (whole cell)  
(ab94173)

**Lane 2 :** 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

#### Secondary

**All lanes :** Goat Anti-mouse IgG (H and L) HRP conjugated at  
1/2500 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors