

p19 INK4d overexpression 293T lysate (whole cell) ab94314

2 图像

概述

| | |
|-------|--|
| 产品名称 | p19 INK4d overexpression 293T裂解物(whole cell) |
| 常规说明 | ab94314 is a 293T cell transfected lysate in which Human p19 INK4d has been transiently over-expressed using a pCMV-p19 INK4d plasmid. The lysate is provided in 1X Sample Buffer. |
| 经测试应用 | 适用于: WB |

性能

| | |
|-----------------|---|
| Mycoplasma free | Yes |
| 形式 | Liquid |
| 存放说明 | Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine) |
| 背景 | Cyclins are important in regulating the cell cycle through their formation of enzymatic complexes with various cyclin dependent kinases. The activity of the cyclin D dependent kinases is controlled by inhibitors such as the INK4 family, which includes INK4a, INK4b, INK4c, and INK4d (p19 INK4d). p19 INK4d is a specific inhibitor of cdk 4 / cdk 6 and is periodically expressed during cell cycle. It is maximally induced as cells enter S phase. |


应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab94314于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|----|------|-------------------------------------|
| WB | | Use at an assay dependent dilution. |

图片



ab94314 at 15µg/lane on an SDS-PAGE gel



All lanes : Anti-p19 INK4d antibody (**ab77223**) at 1/500 dilution

Lane 1 : p19 INK4d overexpression 293T lysate (whole cell)
(ab94314)

Lane 2 : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugated at
1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors