

Pescadillo overexpression 293T lysate (whole cell) ab94281

2 图像

概述

| | |
|-------|--|
| 产品名称 | Pescadillo overexpression 293T裂解物(whole cell) |
| 常规说明 | ab94281 is a 293T cell transfected lysate in which Human Pescadillo has been transiently over-expressed using a pCMV-Pescadillo plasmid. The lysate is provided in 1X Sample Buffer. |
| 经测试应用 | 适用于: WB |

性能

| | |
|-----------------|---|
| Mycoplasma free | Yes |
| 形式 | Liquid |
| 存放说明 | Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine) |
| 背景 | Function: Component of the PeBoW complex, which is required for maturation of 28S and 5.8S ribosomal RNAs and formation of the 60S ribosome. Tissue specificity: Significant levels are detected in a variety of cancer cell lines, including glioblastoma, breast carcinoma, colon carcinoma and cervical carcinoma cells. Levels are abnormally elevated in malignant tumors of astrocytic origin. Similarity: Belongs to the pescadillo family. Contains 1 BRCT domain. PTM: Sumoylated. |


应用

The Abpromise guarantee Abpromise™承诺保证使用ab94281于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|----|------|-------------------------------------|
| WB | | Use at an assay dependent dilution. |

图片



ab94281 at 15µg/lane on an SDS-PAGE gel.



All lanes : Anti-Pescadillo antibody (**ab88543**) at 1/500 dilution

Lane 1 : Pescadillo overexpression 293T lysate (whole cell) (ab94281)

Lane 2 : 293T Non Transfected Lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugated at 1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors