

RFC2 overexpression 293T lysate (whole cell) ab94286

2 图像

概述

| | |
|-------|--|
| 产品名称 | RFC2 overexpression 293T裂解物(whole cell) |
| 常规说明 | ab94286 is a 293T cell transfected lysate in which Human RFC2 has been transiently over-expressed using a pCMV-RFC2 plasmid. The lysate is provided in 1X Sample Buffer. |
| 经测试应用 | 适用于: WB |

性能

| | |
|-----------------|---|
| Mycoplasma free | Yes |
| 形式 | Liquid |
| 存放说明 | Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine) |
| 背景 | Function: The elongation of primed DNA templates by DNA polymerase delta and epsilon requires the action of the accessory proteins proliferating cell nuclear antigen (PCNA) and activator 1. This subunit binds ATP. Disease: Note=RFC2 is located in the Williams-Beuren syndrome (WBS) critical region. WBS results from a hemizygous deletion of several genes on chromosome 7q11.23, thought to arise as a consequence of unequal crossing over between highly homologous low-copy repeat sequences flanking the deleted region. Similarity: Belongs to the activator 1 small subunits family. |


应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab94286于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|----|------|-------------------------------------|
| WB | | Use at an assay dependent dilution. |

图片



ab94286 at 15µg/lane on an SDS-PAGE gel.



All lanes : Anti-RFC2 antibody (**ab88502**) at 1/500 dilution

Lane 1 : RFC2 overexpression 293T lysate (whole cell) (ab94286)

Lane 2 : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugated at 1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors