

Thrombospondin 1 overexpression 293T lysate (whole cell) ab94095

2 图像

概述

产品名称	Thrombospondin 1 overexpression 293T裂解物(whole cell)
常规说明	ab94095 is a 293T cell transfected lysate in which Human Thrombospondin has been transiently over-expressed using a pCMV-Thrombospondin plasmid. The lysate is provided in 1X Sample Buffer.
经测试应用	适用于: WB

性能

Mycoplasma free	Yes
形式	Liquid
存放说明	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: Adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. Binds heparin. May play a role in dentinogenesis and/or maintenance of dentin and dental pulp (By similarity). Ligand for CD36 mediating antiangiogenic properties. Plays a role in ER stress response, via its interaction with the activating transcription factor 6 alpha (ATF6) which produces adaptive ER stress response factors. Similarity: Belongs to the thrombospondin family. Contains 2 EGF-like domains. Contains 1 laminin G-like domain. Contains 1 TSP C-terminal (TSPC) domain. Contains 3 TSP type-1 domains. Contains 8 TSP type-3 repeats. Contains 1 VWFC domain.


应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab94095于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		Use at an assay dependent dilution.

图片



SDS-PAGE - Thrombospondin overexpression 293T lysate (whole cell) (ab94095)

ab94095 at 15µg/lane on an SDS-PAGE gel.



Western blot - Thrombospondin overexpression 293T lysate (whole cell) (ab94095)

All lanes : Anti-Thrombospondin 1 antibody (**ab88529**) at 1/500 dilution

Lane 1 : Thrombospondin 293T Transfected Lysate - Positive Control (ab94095)

Lane 2 : 293T Non Transfected Lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugated at 1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors