

# Mouse Interferon gamma ELISPOT Set ab46586

## 概述

产品名称	小鼠Interferon gamma酶联免疫斑点法(ELISPOT) Set
样品类型	Suspension cells
检测类型	Sandwich (qualitative)
实验步骤	Multiple steps standard assay
种属反应性	与反应: Mouse
产品概述	<p>The Mouse Interferon gamma ELISPOT Set is designed to enumerate cytokine producing cells in a single cell suspension. This method has the advantage of requiring a minimum of in-vitro manipulations allowing cytokine production analysis as close as possible to in-vivo conditions in a highly specific way. This technique is designed to determine the frequency of cytokine producing cells under a given stimulation, and the follow-up of such frequency during a treatment and/or a pathological state. Elispot assay constitutes an ideal tool in the TH1 / TH2 response, vaccine development, viral infection monitoring and treatment, cancerology, infectious diseases, autoimmune diseases and transplantation.</p> <p>This ELISPOT assay is based on sandwich immuno-enzyme technology. Cell secreted cytokines or soluble molecules are captured by coated antibodies avoiding diffusion in supernatant, protease degradation or binding on soluble membrane receptors. After cell removal, the captured cytokines are revealed by tracer antibodies and appropriate conjugates.</p> <p><b>Principal:</b></p> <p>After cell stimulation, locally produced cytokines are captured by a specific monoclonal antibody. After cell lysis, trapped cytokine molecules are revealed by a secondary biotinylated detection antibody, which is in turn recognised by streptavidin conjugated to alkaline phosphatase. PVDF-bottomed-well plates are then incubated with BCIP/NBT substrate. Colored "purple" spots indicate cytokine production by individual cells.</p> <p>Recognizes natural murine Interferon Gamma.</p>
经测试应用	适用于: ELISpot
平台	Reagents

## 性能

## 存放说明

Store at +4°C. Please refer to protocols.

组件	5 x 96 tests	10 x 96 tests	15 x 96 tests	20 x 96 tests
Bovine Serum Albumin	1 x 2g	2 x 2g	3 x 2g	4 x 2g
Clone DB1 Capture antibody	1 x 500µl	2 x 500µl	3 x 500µl	4 x 500µl
Dry Skimmed milk	1 vial	2 vials	3 vials	4 vials
IFN $\gamma$ Biotinylated detection antibody	1 vial	2 vials	3 vials	4 vials
Ready-to-use BCIP/NBT substrate buffer	1 x 50ml	2 x 50ml	3 x 50ml	4 x 50ml
Streptavidin - Alkaline Phosphatase conjugated	1 x 50µl	2 x 50µl	3 x 50µl	4 x 50µl

功能	Produced by lymphocytes activated by specific antigens or mitogens. IFN-gamma, in addition to having antiviral activity, has important immunoregulatory functions. It is a potent activator of macrophages, it has antiproliferative effects on transformed cells and it can potentiate the antiviral and antitumor effects of the type I interferons.
组织特异性	Released primarily from activated T lymphocytes.
疾病相关	In Caucasians, genetic variation in IFNG is associated with the risk of aplastic anemia (AA) [MIM:609135]. AA is a rare disease in which the reduction of the circulating blood cells results from damage to the stem cell pool in bone marrow. In most patients, the stem cell lesion is caused by an autoimmune attack. T-lymphocytes, activated by an endogenous or exogenous, and most often unknown antigenic stimulus, secrete cytokines, including IFN-gamma, which would in turn be able to suppress hematopoiesis.
序列相似性	Belongs to the type II (or gamma) interferon family.
翻译后修饰	Proteolytic processing produces C-terminal heterogeneity, with proteins ending alternatively at Gly-150, Met-157 or Gly-161.
细胞定位	Secreted.

## 应用

**The Abpromise guarantee** **Abpromise™** 承诺保证使用ab46586于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISpot		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors