

Anti-5HT Receptor antibody ab56004

概述

| | |
|-------|---|
| 产品名称 | Anti-5HT Receptor抗体 |
| 描述 | 兔多克隆抗体to 5HT Receptor |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: WB, ELISA |
| 种属反应性 | 与反应: <i>Heliothis virescens</i> |
| 免疫原 | Synthetic peptide corresponding to the C terminal residues of Tobacco budworm [<i>Heliothis virescens</i>] 5-HT receptor (5-hydroxytryptamine receptor). |
| 常规说明 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C. |
| 存储溶液 | <p>pH: 7.2</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 50% Glycerol (glycerin, glycerine), PBS</p> |
| 纯度 | Immunogen affinity purified |
| 纯化说明 | Purified by antigen specific affinity chromatography. |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab56004于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|-------|------|--|
| WB | | Use a concentration of 1 µg/ml. 1µg/ml for 2hrs. |
| ELISA | | Use at an assay dependent concentration. |

靶标

| | |
|------|--|
| 相关性 | The 5-HT receptor (5-hydroxytryptamine receptor) belongs to the G protein coupled receptor 1 family. |
| 细胞定位 | Cell Membrane. Integral membrane protein. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors