

Anti-6X His tag® antibody [3D5] ab206492

重组

1 图像

概述

产品名称	Anti-6X His tag®抗体[3D5]
描述	人单克隆抗体[3D5] to 6X His tag®
宿主	Human
特异性	ab206492 reacts to recombinant proteins containing the 6X His tag® fused to carboxy terminus of fusion proteins. It also has affinity for internal 6X His tag® but no affinity for 6X His tag® fused to amino terminus.
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Species independent
免疫原	Recombinant fragment corresponding to 6X His tag®. (Purified his tagged Fab fragment).
常规说明	This is a recombinant chimeric antibody created for improved compatibility with existing reagents, assays and techniques.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Proclin 300 Constituent: 99% PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	3D5
同种型	IgG1
轻链类型	kappa

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab206492于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.


靶标

相关性 The H-H-H-H-H-H motif is used as a tag on many recombinant proteins to facilitate purification. His-tags can be fused to the amino- or carboxy- termini of proteins in transfected or transformed cells.


细胞定位 Depends upon the localization of the parent protein tagged with hexahistidine.

图片


Why choose a recombinant antibody?




Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-6X His tag® antibody [3D5] (ab206492)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors