

Anti-Ada protein antibody [ADA-1] ab18104

2 References

概述

产品名称	Anti-Ada蛋白抗体[ADA-1]
描述	小鼠单克隆抗体[ADA-1] to Ada蛋白
宿主	Mouse
经测试应用	适用于: RIA, WB, IP
种属反应性	与反应: Escherichia coli
免疫原	Full length protein corresponding to Mouse Ada protein. Please refer to the reference below.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
纯度	Protein A/G purified
克隆	单克隆
克隆编号	ADA-1
同种型	unknown

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab18104于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
RIA		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution. Predicted molecular weight: 40 kDa.
IP		Use at an assay dependent dilution.

靶标

相关性 The methylated ADA protein acts as a positive regulator of its own synthesis, as well as that of other proteins. The transcription-activating function of the ADA protein resides in its N-terminus. It activates the transcription of *alkA*, *alkB* and *aidB*.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors