

Anti-ADCY8/AC8 antibody ab196686

★★★★★ [1 Abreviews](#) [2 References](#) [2 图像](#)

概述

产品名称	Anti-ADCY8/AC8抗体
描述	兔多克隆抗体to ADCY8/AC8
宿主	Rabbit
经测试应用	适用于: IHC-P, ICC/IF
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Synthetic peptide corresponding to Human ADCY8/AC8 (internal sequence). Database link: P40145
阳性对照	Human brain tissue, A549 cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride, 49% PBS PBS is without Mg ²⁺ and Ca ²⁺
纯度	Immunogen affinity purified
纯化说明	Purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen
克隆	多克隆

同种型

IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab196686于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P	★★★★★ (1)	1/50 - 1/1000.
ICC/IF		1/100 - 1/500.

靶标


相关性

This is a membrane-bound, calcium-inhibitable adenylyl cyclase. It may be involved in learning, in memory and in drug dependence.

细胞定位


Cell Membrane; multi-pass membrane protein

图片



Immunofluorescence analysis of A549 cells labeling ADCY8/AC8 antibody with ab196686 at 1/100 dilution. Right panel was preincubated with blocking peptide.

Immunocytochemistry/ Immunofluorescence - Anti-ADCY8/AC8 antibody (ab196686)



Immunohistochemical analysis of paraffin-embedded Human brain tissue labeling ADCY8/AC8 using ab196686 at 1/50 dilution. Right panel was preincubated with blocking peptide.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ADCY8/AC8 antibody (ab196686)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors