

Anti-ADX antibody [EPR4629] ab108257

重组 RabMAb

5 References 2 图像

概述

产品名称	Anti-ADX抗体[EPR4629]
描述	兔单克隆抗体[EPR4629] to ADX
宿主	Rabbit
经测试应用	适用于: WB 不适用于: Flow Cyt or IHC-P
种属反应性	与反应: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
常规说明	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none">- High batch-to-batch consistency and reproducibility- Improved sensitivity and specificity- Long-term security of supply- Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 0.05% BSA, 40% Glycerol (glycerin, glycerine), 59% PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	EPR4629
同种型	IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab108257于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/1000 - 1/10000. Detects a band of approximately 14 kDa (predicted molecular weight: 19 kDa).

应用说明

Is unsuitable for Flow Cyt or IHC-P.

靶标

功能

Participates in the synthesis of thyroid hormones. Transfers electrons from adrenodoxin reductase to the cholesterol side chain cleavage cytochrome P450.


序列相似性

Belongs to the adrenodoxin/putidaredoxin family.
Contains 1 2Fe-2S ferredoxin-type domain.

细胞定位

Mitochondrion matrix.

图片



All lanes : Anti-ADX antibody [EPR4629] (ab108257) at 1/1000 dilution (Purified)

Lane 1 : A549 (Human lung carcinoma epithelial cell) whole cell lysate

Lane 2 : Mouse kidney lysate





Lane 3 : Rat kidney lysate

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution

Predicted band size: 19 kDa

Why choose a recombinant antibody?

 <p>Research with confidence Consistent and reproducible results</p>	 <p>Long-term and scalable supply Recombinant technology</p>
 <p>Success from the first experiment Confirmed specificity</p>	 <p>Ethical standards compliant Animal-free production</p>

Anti-ADX antibody [EPR4629] (ab108257)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors